ResMed

*Changing lives
with every breath*

➤ ResMed Air™ Solutions

⤵ a new beginning...

...in therapy choices.


AirCurve™ 10

CS PACEWAVE

➤ Air Devices: What's new?

Comparison **Air10** series vs **S9** series

➤ Intuitief gebruikersscherm

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> Klinisch menu en schermen

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

Others

> Key features

- Integrated connectivity: telemonitoring

- Oxymetry

Aircurve

➤ Adaptive Servo-Ventilation (ASV)

AirCurve™ 10

CS PACEWAVE

with **Minute Ventilation*** **Adaptive Servo-Ventilation (MV ASV)** therapy

* Minute ventilation = respiratory rate (per minute) x tidal volume of each breath

➤ Indications for use post Serve-HF

- Indications for ASV use have changed as a result of the SERVE-HF study findings. ASV therapy is still indicated to stabilise ventilation in adult patients with CSA, mixed sleep apnoea and periodic breathing, with or without OSA.
- In alignment with regulatory authorities, ASV therapy is now **contraindicated** in patients with chronic, symptomatic heart failure (NYHA 2-4) with reduced left ventricular ejection fraction (**LVEF \leq 45%**) and moderate to severe predominant CSA.


➤ AirCurve 10 CS PaceWave

Overview

The AirCurve 10 CS PaceWave rapidly stabilises breathing because it responds to:

- central apneas, stabilizing ventilation (and therefore blood gases)
and
- obstructive events, stabilizing the upper airway (and therefore maintaining airway patency)

AirCurve 10 CS PaceWave: learns, responds, predicts, optimises

Learning continuously for
personalised therapy

The only ASV therapy to target
the patient's own recent MV

Responding rapidly for
effective therapy

ResMed's most responsive
ASV therapy

Predicting each patient's
unique needs for ease-of-care

Treating your challenging patients
has never been easier

Optimising comfort and
synchrony for compliance

Boost compliance with natural breathing
comfort and advanced synchrony

➤ 1. AirCurve 10 CS PaceWave: learns continuously

To stabilise blood gases it is necessary to stabilise ventilation. To do this AirCurve 10 CS PaceWave **continuously** measures minute ventilation (MV).

- MV target set to 90% of the patient's own ventilation*
- Target MV is continually adjusted to reflect changes in patient's own MV.
- No over-ventilation at any time.


* At the beginning, MV is set to 4l/min and the respiratory rate is set to 15bpm

2. AirCurve 10 CS PaceWave: responds quickly

1. Maintaining MV target

- AirCurve 10 CS PaceWave responds rapidly to changes in minute ventilation
- Monitors recent average minute ventilation (~3-min window)
- Adjusts pressure support up or down as needed to achieve target

3-minute moving window

Increased pressure support when minute ventilation decreases

Decreased pressure support when minute ventilation increases. This prevents the risk of further hyperventilation

2. AirCurve 10 CS PaceWave: responds quickly

1. Stabilizing the upper airway


If the upper airway is closed, ventilation cannot be stabilised.

2 modes:

- In **ASV** mode, fixed EPAP to keep the airway open.

- In **ASV Auto** mode, the algorithm responds to obstructive events


2. AirCurve 10 CS PaceWave: responds quickly

- **Response as and when required:**
 - In **ASV mode**,
 - pressure support (PS) is auto-adjusting and responds instantaneously to changes in minute ventilation
 - EPAP manually set
 - In **ASV Auto mode**,
 - Auto-adjusting PS and
 - Auto-adjusting EPAP that responds on the next breath to flow limitation, snore and obstructive apneas

Response:	Responds to:	Response time:
Pressure support	Minute ventilation	Instantaneously
EPAP	Flow limitation Snore Obstructive apneas	On the next breath

3. AirCurve 10 CS PaceWave: continually predicts

- The AirCurve 10 CS PaceWave algorithm looks at 13 points of each patient's breath, continually mapping the patient's respiratory rate and MV.
- This breath-phase mapping ensures that the AirCurve 10 CS PaceWave provides breaths synchronised with the patient's breathing when the patient's effort decreases.
- When the patient's MV deviates too much from the target, AirCurve 10 CS PaceWave will adjust PS. If target MV is not reached despite PS increase, back-up breaths at a rate of 15 bpm will be provided (if there is no patient RR in memory)


Measurement of 13 points during the breath cycle ensures accurate measurement of minute ventilation and patient-machine synchrony

4. AirCurve 10 CS PaceWave: optimises comfort

RAMP

- EPAP increases linearly over ramp period
 - Start EPAP \Rightarrow Min EPAP
- PS fixed at Min PS during ramp
- Time between 5 - 45 minutes

How does Ramp optimise comfort?

- Heart Failure patients often hyperventilate during treatment onset due to the physical effort exerted before going to bed.
- During Ramp there is no PS reaction to changes in MV and this allows the patient time to relax and lower their minute ventilation.

> The AirCurve 10 CS PaceWave is easy to set-up

Settings	
< Home	
Therapy	
Mode	ASV
EPAP	8.0
Min PS	4.0
Max PS	10.0
Mask	Nasal

3 settings

ASV mode:

- Only 3 therapy settings
- Same default settings provide seamless transition for current ASV users

Settings	
< Home	
Therapy	
Mode	ASVAuto
Min EPAP	4.0
Max EPAP	8.0
Min PS	4.0
Max PS	10.0

4 settings

ASV Auto mode:

- Only 4 therapy settings
- Self-titrating
- Appropriate for healthcare providers wishing to set patients up with auto-adjusting EPAP