Basic 30

- EN Instructions for use
- DE Gebrauchsanweisung
- FR Mode d'emploi
- IT Istruzioni per l'uso
- NL Gebruiksaanwijzing
- SV Bruksanvisning
- FI Käyttöohjeet
- NO Bruksanvisninger
- ES Instrucciones de uso
- PT Instruções de utilização

Technical specifications – Technische Daten – Données techniques –
Dati tecnici – Technische gegevens – Tehniska data – Tekniset tiedot –
Tehniske data – Datos técnicos – Características técnicas

high vacuum

- 85 kPa
- 638 mmHg

high flow

30 l/min.

IPX1

16 kg
35.2 lbs

10 kg
22 lbs

965 x 510 x 480 mm
37.99 x 20.08 x 18.90 inches
290 x 305 x 325 mm
11.42 x 12.01 x 12.8 inches

ISO 9001
ISO 13485
CE (93/42/EEC), IIa

CE

0123

Without jar
Ohne Behälter
Sans récipient
Senza contenitore
Zonder opvangpot

Utan behållare
Ilman astiaa
Uten holder
Sin recipiente
Sem recipiente

106
70
kPa

Operation
Betrieb
Fonctionnement
Funzionamento
Werkung

Drift
Käyttö
Drift
Funcionamiento
Funcionamento

+40
+5
-20
°C

Operation
Betrieb
Fonctionnement
Funzionamento
Werkung

Drift
Käyttö
Drift
Funcionamiento
Funcionamento

75
30
°C

Operation
Betrieb
Fonctionnement
Funzionamento
Werkung

Drift
Käyttö
Drift
Funcionamiento
Funcionamento

+50
-20
°C

Transport / Storage
Transport / Lagerung
Transport / Stockage
Trasporto / Conservazione
Transport / Opslag

Transport / Lagring
Kuljetus / Varastointi
Transport / Lagring
Transporte / Imacenamiento
Transporte / Armazenagem

95
20
°C

Transport / Storage
Transport / Lagerung
Transport / Stockage
Trasporto / Conservazione
Transport / Opslag

Transport / Lagring
Kuljetus / Varastointi
Transport / Lagring
Transporte / Imacenamiento
Transporte / Armazenagem

AC

230 – 240 V 50 Hz 110 W
230 – 240 V 60 Hz 110 W
100 – 120 V 50 Hz 110 W
100 – 120 V 60 Hz 110 W

LATEX FREE

EN	Instructions for use	2
DE	Gebrauchsanweisung	16
FR	Mode d'emploi	30
IT	Istruzioni per l'uso	44
NL	Gebruiksaanwijzing	58
SV	Bruksanvisning	72
FI	Käyttöohjeet	86
NO	Bruksanvisninger	100
ES	Instrucciones de Uso	114
PT	Instruções de utilização	128
A	Technical Documentation (EMC)	142
B	Safety-related checks	146

Caution: Federal US law restricts this device to sale by or on the order of a health professional.

Subject to change

I Congratulations

The Basic 30 is a high-quality secretion aspirator. Its well-proven piston/cylinder system guarantees maximum suction performance and quiet, dependable operation. The simple handling and cleaning as well as the safety features are additional advantages of the Medela secretion aspirators. A comprehensive range of accessories makes the Basic 30 ideally suited to a wide range of medical applications. Contact us – we will be pleased to advise you!

II Intended use

The Basic 30 is indicated for the vacuum extraction, aspiration and removal of surgical fluids, tissues (including bones), gases, bodily fluids or infectious materials from wounds or from the airways of the patient either during surgery or on the ward.

Applications

- General surgery
- Vacuum extraction

For vacuum extraction, the Basic 30 must be used in combination with Medela accessories (suction cups and tubes). Otherwise, Medela can not guarantee the safe functioning of the system.

- Wound drainage

When the Basic 30 is used for wound drainage, it must be ensured that the wound does not dry out and that the dressing is applied tightly. If an unusually high amount of fluids collects in the secretion container (high blood loss), the wound drainage must be interrupted.

- ENT

A replacement device must always be available for patients for whom a breakdown of the device can lead to a critical situation (e.g. patients with dyspnoea or heavy mucous congestion).

Contraindications

Currently no contraindications are known to us.

Important Note

Compliance with proper surgical procedures and techniques is the responsibility of the physician. Each physician must evaluate the appropriateness of the treatment based on his own knowledge and experience.

III Warnings/Safety instructions

The Basic 30 is approved exclusively for the use as described in these instructions. Medela can only guarantee the safe functioning of the system when the Basic 30 is used in combination with the original Medela accessories (collection system, tubes, filter etc.).

The Basic 30 is EMC-tested in conformity with the requirements of IEC 60601-1-2:2007 and can be used in the vicinity of other EMC-tested devices that fulfill the requirements as outlined in the IEC 60601-1-2 standard. Untested HF sources, radio networks or the like can impair the function of the device and should not be operated in combination with the Basic 30.

Please read and observe these warnings and safety instructions before putting into operation.

Warnings

- For use only by medically trained persons who have been adequately trained in suction procedures and in the use of aspirators. Caution: incorrect use can cause pain and injury to the patient.
- Before you plug in the device, please check that your local power supply is the same as the voltage given on the specification plate.
- The device must not be used for suctioning explosive, easily flammable or corrosive liquids
- The connecting tubing supplied with the device must never come into direct contact with the suction area.
- Before cleaning the device, pull the plug out of the mains socket.

Safety instructions

- The Medela Basic 30 is a medical device that requires special safety measures in regard to EMC. It must be installed and put into operation in accordance with the attached EMC information.
- Portable and mobile RF communications equipment can affect medical devices.
- In the case of overflow, inform the internal technical service immediately and perform the tasks in the service manual.
- In each of the following cases, the device must not be used and it must be repaired by Customer Services:
 - if the power cord or the plug are damaged
 - if the device is not functioning perfectly
 - if the device is damaged
 - if the device shows clear safety defects.
- Keep the power cord away from hot surfaces.
- The plug and the ON/OFF switch must not come into contact with moisture.
- Separation from the mains is only assured through the disconnection of the plug and socket connection.
- Never pull the plug out of the mains socket by pulling on the power cord!
- Never leave the device unattended when it is switched on.

These instructions for use must be kept for later reference.

IV Safety tests

For safety tests (STK), it is assumed that the device is serviced and repaired throughout its service life in accordance with the service manual.

Because the Basic 30 is a device of safety class II (EN IEC 60601-1, A1, A2:1995), the safety tests (STK) are confined to visual inspection of the housing and power cord. This test must be carried out before each use.

Devices of protection class II do not have a protective earth conductor, therefore a check of the earth leakage current is not necessary!

The housing of the Basic 30 is made entirely of insulating material. Tests of the housing leakage current using common measuring instruments will therefore not reveal measurable values.

The applied parts of the Basic 30 are connected to the Basic 30 with non-conductive hoses, containers and intermediate hoses.

Even when suctioning a conductive fluid until the overflow protection device activates (9 g/L NaCl; EN ISO 10079-1;1999), measurements of the patient leakage current using common measuring instruments will not reveal measurable values

The Basic 30 does not have patient leakage current, signal input and output parts and functional earth leakage current.

V Occupational health and safety

Guidelines for Medela and agency workshops.

We recommend a hepatitis B vaccination for all persons who have to deal with potentially contaminated devices.

Further safety measures with regard to further vaccinations or other safety precautions are to be observed individually depending on the region.

When devices are returned to workshops for repair, the devices/components have to be considered as being dangerous and potentially infectious.

Such devices must be stored separately, specially marked and at an isolated location until they are disinfected.

VI Important Notes

By the application of medical aspiration devices, vacuum is normally given as the difference (in absolute figures) between absolute pressure and atmospheric pressure or as negative values in kilopascal (kPa).


In this document, the indication of -10 kPa for example always refers to a pressure range in kPa below atmospheric ambient pressure.

VII Guarantee

5 year guarantee after date of delivery for the device. Medela assumes no liability for modified devices, damage or consequential damage caused by improper handling, inappropriate use as well as manipulation by unauthorised persons. Routine checks and service work are only to be carried out by positions authorised by Medela.

VIII Maintenance/Routine check

Routine checks and service work are only to be carried out by positions authorised by Medela. The routine check is to be carried out 1x per year. (see service manual).

X Accessories

Correct and safe functioning can only be guaranteed with Medela accessories.

Fluid Management System

Reusable system		1 L 077.0110
		2 L 077.0120
		3 L 077.0130
		5 L 077.0150
		077.0420
		Size standard for all jars
Disposable system		1.5 L 077.0082
		2.5 L 077.0085
		1.5 L 077.0083
		2.5 L 077.0086
		Without solidifier
		1.5 L 077.0084
		2.5 L 077.0087
		With solidifier


Further accessories

077.0563

077.0571

077.0922

077.0951

077.0762

077.0012

077.0531

077.0723

077.0731

XI Initial Setup

Use only after instruction by qualified personnel.


1

Preparation for use

Check required parts.

- Basic 30
- Tubing set
- Safety device

2


Assembly


- Decision: disposable or reusable system.

For reusable system:

- Mechanical overflow protection must be open. Pull downwards.
- Position jar in clamp holder.
- Connect tubing.


3

Power

- Plug in.


4

Functional check

Ensure that sufficient vacuum is produced.

- Switch device on.

- Set maximum vacuum.

- Seal the end of the patient tubing with your thumb.

- Compare the vacuum.

Max. vacuum according to technical specifications


+ 2,000 m	– 66 kPa	– 495 mmHg
+ 1,000 m	– 75 kPa	– 563 mmHg
+ 500 m	– 80 kPa	– 600 mmHg
0 m	– 85 kPa	– 638 mmHg

OK?

not OK, see Troubleshooting.

OK

5

Use

- Select correct vacuum according to the particular application.

XII Change disposable liner / reusable jar

- on the basis of a visual check

Reusable system

Disposable system

1

Prepare jar.

Prepare liner.


2


- Remove tubing.
- Switch pump off.
- Remove jar.

- Remove patient tubing.
- Plug liner ports.
- Switch pump off.
- Remove liner.

- Unpack new jar, secure in clamp holder.
- Reconnect tubing.


- Unpack new liner, place in jar.
- Reconnect patient tubing.

- Empty jar.

- Dispose of liner in accordance with local guidelines.

3

- Switch on. Vacuum is built up.
- Check vacuum.

XIII Troubleshooting

Motor not running


Check

- Mains voltage
- Switched on
- Fuse

Insufficient vacuum

Check

- Regulator setting
- Tubing
- Plug-in connections
- Overflow protection
- Suction jar and lid (cracks, brittle areas, discolouration)
- Disposable system (cracks, brittle areas, discolouration)

If the fault cannot be rectified, please contact the internal technical department.

XIV Cleaning and sterilisation

After each use, the parts that had been in contact with the aspirated secretions are to be cleaned, disinfected, sterilised or disposed of. Before cleaning the device, pull the plug out of the mains socket.

General notes

- The internal hospital guidelines take precedence.
- Disconnect the Basic 30 from the mains before cleaning/disinfection.
- Clean the device with a moist cloth and non-abrasive agents.
- Suction jar, lid, disposable liner and other accessories are to be treated in accordance with the supplied instruction sheets.
- Wear protective gloves for cleaning/disinfection.
- Dispose of fluids such as blood and secretions and the parts contaminated with them in accordance with internal hospital guidelines.
- Empty the jar or dispose of liner regularly.
- Detailed information according to Medela cleaning instructions 200.2391.

Medela recommends the following surface cleaning agents

Disinfecting agent	Active ingredients	Amount in 100g	Manufacturer
Hexaquart plus (aldehyde-free) [concentrate]	Didecyltrimethylammonium chloride Lauryl dipropylene triamine Biguanidine acetate	6 g 5.5 g 4.2 g	B. Braun Medical AG www.bbraun.ch
Perform [concentrate]	Potassium peroxomonosulphate Sodium benzoate Tartaric acid Non-anionic surfactants Soap phosphate Perfumes	20 g 15 g 10 g 5 – 15% 5 – 15% < 5%	Schülke & Mayr www.schuelkemayr.com

	Cleaning	Disinfection	Sterilisation
Pump housing, Mains cable	✓ Wipe off with a damp cloth.	✓ With disinfecting agent (see above)	✗  Can not be autoclaved
Disposable products (e.g. disposable liner)	 These are disposable items and must not be reused. Replace after every use! Reprocessing could cause loss of mechanical, chemical and/or biological characteristics.		

Cleaning the reusable jar and lid (PSU and PC)

For optimum user safety, we recommend checking the reusable jar and lid before each use as natural damage can occur with time. This check consists of a visual check and evacuation of the jar at maximum vacuum. Faulty jars and lids should be replaced.

Polysulfone (PSU) jar

Disinfect and sterilise for 10 minutes at 134°C/273°F.

Remove all residues from chemical disinfection agents before autoclaving. Do not stack when autoclaving. Do not use any phenol-based cleaning agents. With proper use, reusable jars (PSU) and lids have a service life of 2 years or 200 autoclaving cycles.

After this time, jars and lids must be replaced for safety reasons.

Polycarbonate (PC) jars

Do not autoclave. The maximum disinfection temperature is 100°C/212°F. Information from the manufacturer regarding concentration, soaking times or temperature must be observed.

Silicone tubing

Can be autoclaved for 10 minutes at 134°C/273°F.

XV Disposal

The Basic 30 comprises metals and plastics and should be disposed of in accordance with the European directives 2002/95/EC and 2002/96/EC. Additional, local guidelines must also be observed.

XVI Article numbers and descriptions

Correct and safe functioning can only be guaranteed with Medela accessories.
Optional accessories - not included in the scope of delivery

Reusable system

077.0110 1 L PSU jar

077.0120 2 L PSU jar

077.0130 3 L PSU jar

077.0150 5 L PSU jar

077.0420 Lid large, conical patient connection

077.0125 0.25 L PSU jar

077.0155 0.5 L PSU jar

077.0440 Lid small, conical patient connection

Disposable system

077.0082 1.5 L PC jar

077.0085 2.5 L PC jar

077.0083 1.5 L disposable liner

077.0086 2.5 L disposable liner

077.0084 1.5 L disposable liner with solidifier

077.0087 2.5 L disposable liner with solidifier

Further accessories

077.0563 Change-over-valve

077.0571 Overflow protection/
bacteria filter

077.0922 Silicone tubing
Ø 7 x 12 mm, 60 cm
with 2 coupling pieces

077.0951 Disposable tubing
(PVC), 180 cm with
fingertip, sterile

077.0762 Suction cup SilcCup,
Ø 60 mm

017.0012 Apgar Timer

077.0531 Quiver with holder

077.0723 Foot On/Off switch

077.0731 Foot vacuum regulator

A

Accessories 2, 3, 7, 11, 13, 14
Ambient pressure 5
Applications 2
Article numbers 13

C

Check 3, 4, 5, 9, 10, 12
Cleaning 2, 3, 11, 12
Consequential damage 5
Contraindications 2

D

Disinfecting agent 11
Disinfection 11, 12
Disposable system 7, 8, 10, 11, 13
Disposal 12

E

ENT 2

F

Foot switch 6
Functional check 9
Functioning 2, 3, 7, 13
Fuse 11

G

Guarantee 5
Guidelines 4, 10, 11, 12

H

Handling 2, 5

I

Indication 5
Initial setup 8

J

Jar 8, 10, 11, 12, 13

K

Kilopascal 5

L

Liability 5
Liner 10, 11, 12, 13

M

Mains 3, 11
Maintenance 5
Mobile stand 6

O

ON/OFF switch 3, 6
Overflow protection 11, 14

P

Patient tubing 9, 10
Plug 3, 8, 11
Pressure 5
Putting into operation 3

R

Range of accessories 2
Reusable jar 6, 10
Reusable system 7, 8, 10, 13
Routine check 5

S

Safety features 2
Safety instructions 3
Safety measures 4
Safety set 6, 8
Scope of delivery 13
Silicone tube 6
Silicone tubing 12, 14
Standard rail 6
Sterilisation 11, 12
Suction area 3

T

Troubleshooting 9, 11
Tubes 2, 3
Tubing 3, 8, 9, 10, 12, 14, 34
Tubing set 8

U

Use 2, 3, 5

V

Vacuum 2, 5, 9, 10, 12, 14
Vacuum gauge 6
Vacuum regulator 6, 14
Voltage 3, 11

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Electromagnetic Emissions

Basic 30 is intended for use in the electromagnetic environment specified below. The customer or the user of Basic 30 should assure that they are used in such environment.

Emission Tests	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	Basic 30 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference with nearby electronic equipment.
RF emissions CISPR 11	Class B	Basic 30 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 60000-3-3	Complies	

Warning – Basic 30 should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, Basic 30 should be observed to verify normal operation in the configuration in which it will be used.

IEC 60601-1-2, Table 202

Electromagnetic Immunity

Basic 30 is intended for use in the electromagnetic environment specified below. The customer or the user of Basic 30 should assure that they are used in such environment.

Immunity Tests	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input / output lines	± 2 kV for power supply lines ± 1 kV for input / output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % UT (>95 % dip in UT) for 0.5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 sec	<5 % UT (>95 % dip in UT) for 0.5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of Basic 30 requires continued operation during power mains interruptions, it is recommended that Basic 30 be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels of a typical commercial or hospital environment.

Electromagnetic Immunity

Basic 30 is intended for use in the electromagnetic environment specified below. The customer or the user of Basic 30 should assure that Basic 30 is used in such environment.

Immunity Tests	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of Basic 30, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m) Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range.^b Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4	3 V/m 80 MHz to 2.5 GHz	3 V/m	

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed RF transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which Basic 30 are used exceeds the applicable RF compliance level above, Basic 30 should be observed to verify normal operation. If abnormal operation is observed, additional measures may be necessary, such as reorienting or relocating Basic 30.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

IEC 60601-1-2, Table 206

Recommended separation distance between portable and mobile RF communications equipment and Basic 30			
Basic 30 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of Basic 30 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and Basic 30 as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter	Separation distance according to frequency of transmitter		
W	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

B Safety Tests

DIN EN 62353: 2008-8

Hospital	Basic 30-Serial number
Location	

Visual checks	NO	YES
Cracks, brittle areas, discoloration on housing and display		
Basic 30 is clean with no dirty areas		

Measurements*	OK	Not OK
Max Vacuum (use Basic 30 with 1 L jar and tubing set) Set vacuum to - 85 kPa (-638 mmHg). Using a short piece of suitable tubing (7 mm ID) connect suction tube to vacuum gauge and start unit, wait 1 minute for the pump to reach equilibrium	☐ Measured value is between -85 kPa (-638 mmHg) +/- 5%	☐ Measured value _____
Max. Flow (keep settings from above) Disconnect vacuum gauge and connect flow meter instead, wait 2 minutes for the pump to reach equilibrium	☐ ~ 30 L/min	☐ Measured value _____

Overall evaluation	Faults	Measures
No safety or functional technical faults were detected.	☐	Basic 30 is ready for use
Faults that are not safety relevant were detected. No direct risk.	☐ 1 _____ 2 _____ 3 _____	Basic 30 is ready for use, Faults must be corrected by _____
Safety relevant faults were detected.	☐ 1 _____ 2 _____ 3 _____	Do not use Basic 30 until the faults are corrected. Send Basic 30 to Medela Service center.

* Required Test Equipment
Vacuum gauge Class 2.5 % or better (074.0006)
Flow Meter air 20°, 1 bar abs, Class 4 % or better

Next test

Tester/Signature	Place/Date
------------------	------------

Please see accompanying papers
 Begleitpapiere beachten
 Lire les documents joints à l'appareil
 Attenersi alla documentazione
 Bigjevoegde documenten in acht nemen
 Palun vaadake kaasas olevaid dokumente
 Katso mukana toimitetuista asiakirjoista
 Vennligst se vedlagte dokumenter
 Aténgase a la documentación adjunta
 Verifique os documentos de transporte

IPX1

Splash-proof
 Tropfwassergeschützt
 Étanche à la pluie
 Protezione dall'acqua
 Spatwaterdicht
 Stänksäker
 Kosteaa tila
 Dryppssikker
 Antigoteo
 À prova de salpicos

Protection Class II
 Schutzklasse II
 Classe de protection II
 Classe di protezione II
 Beschermingsklasse II
 Skyddsklass II
 Suojaluokka II
 Beskyttelse klasse II
 Clase de protección II
 Classe de protecção II

Potential equalization
 Potentialausgleich
 Equilibrage de potentiel
 Compensazione del potenziale
 Potentiaalvereffening
 Potentialutjämning
 Potentiaalintasaus
 Spenningsutjevning
 Ecualización potencial
 Igualização de tensão

AP class appliance
 Geräte der Klasse AP
 Appareil de classe AP
 Apparecchio di classe AP
 Klasse AP apparaat
 Apparat, klass AP
 AP-luokan laite
 Apparat i AP-klassen
 Aparato de la clase AP
 Aparelho de classe AP

Type CF
 Typ CF
 Type CF
 Tipo CF
 Type CF
 Type CF
 Tyyppiä CF
 Type CF
 Tipo CF
 Tipo CF

Autoclavable e.g 134°C/10Min.
 Autoklavierbar z.B. 134°C/10Min.
 Autoclavable p. ex. 134°C/10Min.
 Autoclave, p.es., 134°C/10Min.
 Autoclaveerbaar, bijv. 134 °C/10min.
 Går att autoklavera, t.ex. 134°C/10min.
 Autoklaavattava esim. 134 °C/10min. lämpötilassa
 Kan autoklaveres ved 134°C/10min.
 Esterilizable en autoclave p. ej. 134° C/10min.
 Autoclavável p. ex., a 134 °C/10Min.

Not autoclavable
 Nicht autoklavierbar
 Non autoclavable
 Non in autoclave
 Niet autoclaveerbaar
 Går inte att autoklavera
 Ei autoklaavattava
 Kan ikke autoklaveres
 No esterilizable en autoclave
 Não autoclavável

Disposable article – do not reuse
 Einweg-Artikel – nicht wiederverwendbar
 Article à usage unique – ne pas réutiliser
 Articolo monouso – non riutilizzare
 Wegwerpartikel – niet opnieuw bruikbaar
 Engångsartikel – får ej återanvändas
 Kertäkäyttöartikkeli – ei uudelleen käytettävä
 Engangsprodukt – ikke bruk på nytt
 Artículo desechable – no reutilizar
 Artigos descartáveis – não reutilizar

International Sales

Medela AG, Medical Technology
Lättichstrasse 4b
6341 Baar
Switzerland
Phone +41 (0)41 769 51 51
Fax +41 (0)41 769 51 00
ism@medela.ch
www.medela.com

Australia

Medela Pty Ltd,
Medical Technology
3 Arco Lane, Heatherton
Vic 3202
Australia
Phone +61 (0) 3 9552 8600
Fax +61 (0) 3 9552 8699
contact@medela.com.au
www.medela.com.au

Canada

Medela Inc.
4160 Sladeview Crescent Unit # 8
Mississauga, Ontario
Canada, L5L 0A1
Phone +1 905 608 7272
Fax +1 905 608 8720
info@medela.ca
www.medela.ca

China

Medela (Beijing)
Medical Technology Co., Ltd.
Unit M, 4F,
No 1068 Wuzhong Road,
Minhang District
Shanghai 201103
China
Phone +86 21 62368102
Fax +86 21 62368103
info@medela.cn
www.medela.cn

France

Medela France Sarl
14, rue de la Butte Cordière
91154 Etampes cedex
France
Phone +33 (0)1 69 16 10 30
Fax +33 (0)1 69 16 10 32
info@medela.fr
www.medela.fr

Germany

Medela Medizintechnik
GmbH & Co. Handels KG
Postfach 1148
85378 Eching
Germany
Phone +49 (0)89 31 97 59-0
Fax +49 (0)89 31 97 59 99
info@medela.de
www.medela.de

Italy

Medela Italia Srl
Via Turrini, 19 - Loc. Bargellino
40012 Calderara di Reno (BO)
Italy
Phone +39 051 72 76 88
Fax +39 051 72 76 89
info@medela.it
www.medela.it

Japan

Medela K.K.
KDX Kiba Building 3F
5-12-8, Kiba, Koto-ku,
Tokyo 135-0042
Japan
Phone +81 3 3820-5187
Fax +81 3 3820-5181
info@medela.jp
www.medela.jp

Netherlands & Belgium

Medela Benelux BV
Uilenwaard 31
5236 WB 's-Hertogenbosch
Netherlands
Phone +31 73 690 40 40
Fax +31 73 690 40 44
info@medela.nl
info@medela.be
www.medela.nl
www.medela.be

Poland

Medela Polska Sp. z o.o.
Ul. Lewinowska 8
03-684 Warszawa
Poland
Phone +48 22 864 38 65
Fax +48 22 865 12 50
biuro@medela.pl
www.medela.pl

Russia

OOO Medela
Vavilova str. 97, entrance 3
117335 Moscow
Russia
Phone +7 (495) 980 6194
Fax +7 (495) 980 6194
info@medela-russia.ru
www.medela-russia.ru

Spain & Portugal

Productos Medicinales
Medela, S.L.
c/ Manuel Fernández Márquez, 49
08918 Badalona (Barcelona)
Spain
Phone +34 93 320 59 69
Fax +34 93 320 55 31
info@medela.es, info@medela.pt
www.medela.es, www.medela.pt

Sweden

Medela Medical AB
Box 7266
187 14 Täby
Sweden
Phone +46 (0)8 630 09 40
Fax +46 (0)8 630 09 48
info@medela.se
www.medela.se

Switzerland

Medela AG, Medical Technology
Lättichstrasse 4b
6341 Baar
Switzerland
Phone +41 (0)848 633 352
Fax +41 (0)41 769 51 00
contact@medela.ch
www.medela.ch

United Kingdom

Medela UK Ltd.
Huntsman Drive
Northbank Industrial Park
Irlam, Manchester M44 5EG
UK
Phone +44 870 950 5994
Fax +44 870 389 2233
info@medela.co.uk
www.medela.co.uk

USA

Medela Inc.
P.O. Box 660
1101 Corporate Drive
McHenry, IL 60050
USA
Phone +1 877 735 1626
Fax +1 815 363 2487
suction@medela.com
www.medelasuction.com